VISCOSIMETROS ROTACIONALES ROTACIONAL VISCOSIMETERS

ST-2020-L 1001616
ST-2020-R 1001617

Indice

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Index

1. Safety	12
2. Packaging List	12
3. Introduction	12
3.1 Measurement method	12
4. Technical Specifications	13
4.1 Measurement range for every model	13
5. Installation	14
5.1 Viscosimeter assembly	14
5.2 Power supply connection.....	14
6. Operation.....	15
6.1 Control panel description.....	15
6.2 Information shown in the display	15
6.3 Spindle assembly.....	15
6.4 Temperature probe assembly	16
6.5 Start measurement configuration	16
7. Instructions for measurement	17
7.1 Sample preparation.....	17
7.2 Select the correct Full Scale.....	17
7.3 Spindle/Speed selection table	18
7.4 Details for a better spindle selection	19
8. Printer connection and RS-232 communication	20
9. General precautions and advices	20
Adapter for small sample volume 1000996	22
Adapter assembly:.....	22

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1. Safety

Please bear in mind the safety information of the liquids used.

See also chapter 9.

2. Packaging List

- Briefcase (Two cases on the 1001617)
- Main unit
- Rack shaft
- Support feet.
- Spindle Model «L» # 1; # 2; # 3, # 4
(Placed on the same case)
Model «R» R2; R3; R4; R5; R6; R7
(Placed on a separate case)
- Spindle binder
- Spindle protector
- Spanner
- Power Inlet: 110-230 VAC Outlet: 15V 1.5A
- Temperature probe
- Cable RS-232

3. Introduction

The rotational viscosimeter is used for measuring dynamic viscosity. This feature expresses the substances resistance to distortion when a cutting effort is being applied.

This is a highly dependent property of the measurement temperature.

As an example, we have the chocolate: in a small temperature margin, its almost infinite value viscosity (at about 10°C is completely solid) at a very low value (at about 60°C is completely liquid).

Viscosity is measured by establishing the necessary force (cutting or shear effort (N/m²) to move the particles of material g with a particular distortion speed (speed gradient, (S-1)). The viscosity (h) is obtained as a result of the ratio between the cutting effort and the speed gradient.

The measurement unit is the Poise (P), although in practice we use either its hundredth part (cP) or mPa.s in the International System.

The relationship with static viscosity is as follows:

Dynamic V = Static V x density

Other concepts to bear in mind on dynamic viscosity measurement are:

Thixotropy: variation of viscosity with time. For instance: «Quick dry» glues and paints.

Newtonian / Non-Newtonian: constant viscosity in time / varying viscosity

3.1 Measurement method

Figure 3.1 offers a schematic of the measurement system.

The system is made up of a motor (1) which makes the spindle (2) turn, immersed in the sample (3), and joined by a spring (4).

Due to the resistance the sample offers to the spindle, its turn is retarded with regard to the motor.

Rotational viscosimeters use this «retard» to calculate dynamic viscosity.

This «retard» increases the spindle size and the angular speed (rpm).



Fig. 3.1

4. Technical Specifications

<u>Power supply:</u>	15V DC 1.2A
<u>AC / DC supply:</u>	Inlet 100-240 VAC 50/60Hz
<u>Dimensions main unit: height/width/depth</u>	190 x 100 x 145 mm
<u>Dimensions assembled:</u>	620 x 360 x 230 mm
<u>Spindle number:</u>	Model «L» # 1; # 2; # 3, # 4
	Model «R» R2; R3; R4; R5; R6; R7
	R1 is sold apart.
<u>Error:</u>	± 2%
<u>Repeatability:</u>	0.5%

4.1 Measurement range for every model


Suitable for ST2020-R

ST-2020- L

Measurement range (mPa.s) = (cP) 1 ... 600.000

Speed range (rpm) 1 ... 60 by increases of 1rpm


ST-2020- R

Suitable for ST2020-L


Measurement range (mPa.s) = (cP) 1 ... 40.000.000

Speed range (rpm) 0.1 ... 99.9 by increases of 0.1rpm

5. Installation

5.1 Viscosimeter assembly

«Yellow cap» for motor shaft protection.

Keep always the motor shaft protected.

ATTENTION

Please do all the assembly operations with the «Yellow cap» on to protect the motor shaft.

During the main unit manipulation, keep always the viscosimeter shaft protected.

5.2 Power supply connection

The viscosimeter does not connect directly to the mains. A charger with a 15DV 1.5A outlet must be used.

The charger must be connected to the mains by means of a switch without earth wire.

Use only the charger supplied by J.P. SELECTA, s.a.

You can use a similar one but be careful with the connector polarity.

Connector polarity:

Exterior: negative.

Inner: positive.

6. Operation

6.1 Control panel description



Fig 6.3.1.

Resets all parameters to the default value.

«Increase» Key

Prints the display result.
See 6.8

Moves the display cursor to the left.

Begins the measurement with the indicated parameters.

Selects the spindle.
Every pulse changes to no1, no2...

«Decrease» Key

Speed selection


Fig 6.3.2.

6.2 Information shown in the display

Probe measure temperature

Selected speed(rpm)

% Full Scale.
See 6.4

Selected spindle

Viscosity measured

Measurement units

20.0°C 60.0RPM 45.7%
SP 4
8888888.8 mPa.s

6.3 Spindle assembly

The spindle assembly is an extremely careful operation.

An incorrect assembly could damage the motor shaft.

Please, follow the next instructions thoroughly:

a) The spindle protector must be assembled.

To assemble the spindle protector: (See fig 6.3.1)

1- Insert both protector holes in the viscosimeter supports.

2- Turn the protector clockwise.

b) Remove the «yellow cap» (Motor shaft protector).

c) Take the motor shaft and push upwards.

d) By keeping the shaft upwards, fit the spindle and then turn counter-clockwise.
(See fig 6.3.2)

Set spindle:

1- Push the shaft upwards

2- Fit the spindle

3- Turn counter-clockwise


ATTENTION

PLEASE BE CAREFUL WHEN FITTING THE SPINDLE.
THE MOTOR SHAFT IS A DELICATE PART AND DIFFICULT TO REPLACE.

ATTENTION

PLEASE AVOID HITTING ON THE SPINDLE SIDE.
ALWAYS USE THE SPINDLE PROTECTOR.

Temperature probe connector

**6.4 Temperature probe assembly**

It is done once the spindle is assembled.

An incorrect assembly could definitely damage the motor shaft.

Please follow carefully the next instructions:

- First, the spindle protector must be assembled.

To assemble the spindle protector: (See fig 6.3.1)

- Insert both protector holes in the viscosimeter supports.
- Turn the protector clockwise.

- Remove the «yellow cap» (Motor shaft protector).
- Take the motor shaft and push upwards.
- By keeping the shaft upwards, fit the spindle and then turn counter-clockwise. (See fig 6.3.2)
- Release the spindle.

6.5 Start measurement configuration

Once the spindle and the temperature probe are fitted:

Switch on the main switch. (On the rear part of the main unit).

2. Push «RUN» to begin measurement.

3. Push «RESET» to stop measurement.

Speed selection:

1. Push «SPEED».

2. Insert units and tens for speed by using «◀» or «▶» keys.

3. Select «▲» or «▼» keys.

4. Push «SPEED» to validate the entry and then exit.

Measure spindle selection:

1. Push «SPINDLE» key to select the spindle number which has been assembled.

2. Push «SPINDLE» successively till the display shows the spindle number used.

It's not necessary to validate any key.

7. Instructions for measurement

To correctly measure the viscosity, you must use the suitable material, control the ambient temperature and select the correct measurement scale.

7.1 Sample preparation

The liquid must be placed in a 600ml «Low-form» beaker.

You need 500ml of sample.

Place the beaker under the viscosimeter and move the main unit down till the spindle mark arrives to the liquid level. See Fig. 7.1


Fig 7.1 Adjust the viscosimeter height till the liquid level arrives to the spindle mark.

7.2 Select the correct Full Scale

Any combination spindle / speed (rpm) establishes a viscosimeter measurement range.

Spindle effect on Full Scale:

Big spindles, number 1 & 2, are used to measure low viscosity values.

Small spindles, number 3 & 4, are used to measure high viscosity values.

Speed effect (rpm) on Full Scale (FS):

High speed is used to measure low viscosity values.

Low speed is used to measure high viscosity values.

Bearing in mind these effects, any combination speed / spindle has a different Full Scale (FS).

When selecting a speed and a spindle, the display shows the Full Scale.

For instance: (See table 7.2.1)

Combination: spindle = no. 1; 60 rpm has a FS = 100cP

and combination: spindle = no. 1; 6 rpm has a FS = 1000cP

If the sample viscosity value is absolutely unknown:

- Put spindle no. 2
- Begin with a speed of 10 rpm.
- If the display indicates EEEEEE, it means that the sample has a higher viscosity than the FS. Change then to a smaller spindle. (No. 3 or no. 4)
- Check the value on the display upper right corner. (See figure 7.2) This value indicated the measurement viscosity percentage on the Full Scale that corresponds to the combination spindle/rpm.
- This percentage should be as closest as possible to 100%.
- Modify speed or the spindle as appropriate.

Indication % Full Scale

Fig 7.2 Percentage Full Scale.

See chapter 7.3 concerning other details for spindle and speed selection.

7.3 Spindle/Speed selection table

These tables are of great help for selecting the suitable combination spindle/speed. Only several speeds are displayed.

The equipment display shows the FS value during speed or spindle selection.

Table 7.3.1 Full Scale: Model ST-2020 L (cPs)

rpm	Spindle			
	1	2	3	4
1	6.000	30.000	120.000	600.000
6	1.000	5.000	20.000	100.000
10	600	3.000	12.000	60.000
20	300	1.500	6.000	30.000
30	200	1.000	4.000	20.000
40	150	750	3.000	15.000
50	120	600	2.400	12.000
60	100	500	2.000	10.000

Tabla 7.3.2 Full Scale: Model ST-2020 R (cPs)

rpm	Spindle						
	R1	R2	R3	R4	R5	R6	R7
0.1	100.000	400.000	1.000.000	2.000.000	4.000.000	10.000.000	40.000.000
0.3	33.333	133.333	333.333	666.667	1.333.333	3.333.333	13.333.333
0.5	20.000	80.000	200.000	400.000	800.000	2.000.000	8.000.000
0.6	16.667	66.667	166.667	333.333	666.667	1.666.667	6.666.667
1.0	10.000	40.000	100.000	200.000	400.000	1.000.000	4.000.000
1.5	6.667	26.667	66.667	133.333	266.667	666.667	2.666.667
2.0	5.000	20.000	50.000	100.000	200.000	500.000	2.000.000
2.5	4.000	16.000	40.000	80.000	160.000	400.000	1.600.000
3.0	3.333	13.333	33.333	66.667	133.333	333.333	1.333.333
4.0	2.500	10.000	25.000	50.000	100.000	250.000	1.000.000
5.0	2.000	8.000	20.000	40.000	80.000	200.000	800.000
6.0	1.667	6.667	16.667	33.333	66.667	166.667	666.667
10.0	1.000	4.000	10.000	20.000	40.000	100.000	400.000
12.0	833	3.333	8.333	16.667	33.333	83.333	333.333
20.0	500	2.000	5.000	10.000	20.000	50.000	200.000
30.0	333	1.333	3.333	6.667	13.333	33.333	133.333
50.0	200	800	2.000	4.000	8.000	20.000	80.000
60.0	167	667	1.667	3.333	6.667	16.667	66.667
100.0	100	400	1.000	2.000	4.000	10.000	40.000

7.4 Details for a better spindle selection

There are more criterions to follow to select the appropriate combination spindle/speed, besides what we have explained on chapter 7.2.

There are two starting criterions:

Measure it with the slightest measurement error possible.

Measure to be able to compare with other measures.

7.4.1 Measure it with the slightest measurement error possible

We have to find a suitable combination spindle/speed that comply the following requirements (in order of importance):

1. The average value (shown in the display) would be as close as possible to the Full Scale. That is, that the percentage will be the highest possible.

2. The speed would be the highest possible.

For instance:

- We work with ST-2020-L, and apply table 7.3.1

- We want to measure a sample of about 5000cP.

- From table 7.3.1, we could use one of the following combinations:

- Spindle no 1; speed 1 rpm

- Spindle no 2; speed 5 rpm (it does not appear on the table)

- Spindle no 3; speed 20 rpm

- Spindle no 4; speed 50 rpm

- Spindle no 4; speed 60 rpm

and others.

In this case, we will use the next combination:

Spindle no 3; speed 20 rpm

As it is the higher speed (rpm) and higher FS percentage.

7.4.2 Measure to be able to compare among other equipments

In many occasions (for instance, in quality control), the measurements made must be compared to the ones made in other laboratories.

When indicating a measurement result, we must indicate all the conditions such as: viscosity, temperature and spindle number.

«Newtonian» and «Non-newtonian» fluids:

There are two kinds of fluids:

- «Newtonian»: The dynamic viscosity measurement is **not** affected when changing the average speed or the spindle.

- «Non-newtonian»: The dynamic viscosity measurement **do** affect when changing the average speed or the spindle.

For instance: A «Newtonian» liquid of 5000cP provides the same reading in the viscosimeter if measured under the following two conditions:

- Spindle no 3; speed 20 rpm

- Spindle no 4; speed 50 rpm

If the liquid would be «Non-newtonian», the reading would be different.

In practice, most of the liquids are «Non-newtonian». Only some silicone oils are «Newtonian». That's why you have to enter the spindle number and the speed used.

8. Printer connection and RS-232 communication

You can connect a serial printer as 4120113 to print the data.

You can also ask J.P. Selecta for the software to directly download the measurements to the PC.

The serial communication follow the next parameters:

9600bps, 8 data bits, 1 stop bit, non-parity.


Cable schematics:

9 pin male	9 pin male
2 -----	3
3 -----	2
5 -----	5

9. General precautions and advices

- Keep the motor shaft always protected with the yellow cap.
- Assemble and dismantle the spindle carefully. The motor shaft is a very delicate part.
- It's not necessary to «hold tight» the spindle to the motor shaft.
- Do not use any tool to «tighten» the spindle to the motor shaft.
- Do not turn the spindle manually once installed in the viscosimeter.
- Keep the spindle and the motor shaft threads clean.
- Once the spindle is assembled, do not lie the viscosimeter down.
- To turn the knob upwards and, hold the viscosimeter with the other hand.
- When changing the spindle, remember to select it on the viscosimeter.
- When changing the fluid, you must clean the spindle and its protector.
- Do not use other viscosimeters spindles.
- Ambient temperature affects the viscosimeter measurement a lot.

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



Adapter for small sample volume 1000996

When we have only a small sample volume, we must use the adapter 1000996.

The Sample Adapter consists of:

- Extension link for spindle with coupling nut to connect to the motor axle.
- Special spindles (spindle no. 21, 27, 28 and 29)
- Sample chamber with thermostating jacket.
- Adjustable base.

Adapter assembly:

1. Fix the extension link to the motor axle. **ATTENTION: See 6.3 concerning how to do this action in order not to damage the motor.**

2. Hang the selected spindle: no. 21, no. 27, no. 28 or no. 29.

Lift the viscosimeter up as much as you can in order to facilitate the operation.

3. Place the chamber in the centre of the adapter base and fix it with three screws.

Do not excessively overtighten it as the position must be corrected later on to centre the adapter to the viscosimeter shaft.

Lower down the viscosimeter slowly by turning the knob.


Move the adapter base so that the spindle can smoothly enter into the sample chamber.

Adjust the adapter base level with the help of the three screws and the level bubble.

Fill the chamber with 10 - 20ml liquid sample (this operation can be made before or after inserting the spindle).

If you want to temper the sample, connect a recirculation temperature bath to the in and out water jackets on the adapter side.

4. Configure, in the viscosimeter panel, the spindle number. Select the viscosity and start the measurement by pressing the «Test» key (reading must be between 20 and 90%).

Adapter for small sample volume 1000996

When we have only a small sample volume, we must use the adapter 1000996.

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

**Notice to customers:**

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.