


## 組織均質頭 Dispersing Elements



Model	Capacity (ml)	Dimension (mm)	Application
K-7S	10	Ø7 x 130	微量細胞組織之均質、乳化、混合 Applies to homogenize and mix small volumes of cell tissue
K-12S	250	Ø12 x 170	動、植物組織之破碎 懸浮物及乳劑之處理 乳化、混合、萃取等作用 Applies to disperse plant and animal tissue Applies to dispose suspensions and emulsions With homogenizing, mixing and extracting actions
K-20S	2000	Ø20 x 200	動、植物組織之破碎 懸浮物及乳劑之處理 乳化、混合、萃取等作用 Applies to disperse plant and animal tissue Applies to dispose suspensions and emulsions With homogenizing, mixing and extracting actions
K-30	3000	Ø30 x 250	有機與無機物之破碎 乳化、混合、萃取等作用 Applies to disperse organic and inorganic matter With homogenizing, mixing and extracting actions